NEW ADVANCEMENT IN BOVINE REPRODUCTION

VetMott, Inc.

MULTI (850µL, 3ML) SPERM SEPARATION DEVICES

VetMotl Sperm Separation Devices, first-of-their-kind lab tools for use in veterinary assisted reproductive technology procedures.


Our devices have been commercialized for use in human fertility clinics as ZyMōt® Sperm Separation devices since 2015.



VMB0850 (850µl) and VMB03000 (3ml)



blastocyst rate



VetMotl prepared spermatozoa induced a significantly higher blastocyst rate at day 7 compared to DGC¹.

- ✓ Optimized for use in bovine IVF
- ✓ Conventional or sexed semen
- No special media needed
- Replaces traditional methods including swim up, wash, and density gradient centrifugation (DGC)
- ✓ Ready to use & easy to adopt

FOR MORE INFORMATION:

Kristopher Scott kscott@vetmotl.com • 678-836-8431

www.vetmotl.com





401 Professional Drive Suite 130 Gaithersburg, MD. 20879-3429 240-454–9893

Sources: 1Lorenzen 2022